



## XRF Platinum Program for Leco



### XRF Crucibles suitable for Leco

Diameter internal		Height (H)	Base thickness (s)	Weight approx.	Heraeus drawing
top (OD)	bottom (UD)				
mm	mm	mm	mm	g	
37	20	35	0,5	36	80623536



### Dispersion Hardened (DPH) Alloy


Through a special process, Heraeus has created a Dispersion Hardened (DPH) alloy.

The inclusion of finely dispersed zirconia offers an economical way to increase the mechanical strength of the product resulting in longer service times.

### Engraving possibilities




For easier tracking of your platinum products in the lab Heraeus offers different engraving possibilities.



### XRF Casting dishes suitable for Leco (round version)

Diameter internal		Flange diameter (FD)	Height (H)	Base thickness (s)	Weight approx.	Heraeus drawing
top (OD)	bottom (UD)					
mm	mm	mm	mm	mm	g	
32,0	30,0	39	6,0	0,8	21	80623259
32,0	30,0	39	6,0	1,0	27	80623259
32,6	31,5	41	3,5	1,0	29	80668768
34,0	32,0	41	6,0	0,8	23	80623267
34,0	32,0	41	6,0	1,0	29	80623267
34,0	32,0	41	6,0	1,2	36	80623267
35,6	34,5	44	3,5	1,0	36	80668782
37,0	35,0	44	6,0	0,8	26	80623268
37,0	35,0	44	6,0	1,0	33	80623268
37,0	35,0	44	6,0	1,2	40	80623268
40,0	38,0	48	6,0	0,8	34	80623270
40,0	38,0	48	6,0	1,0	43	80623270
40,6	39,5	49	3,5	1,0	44	80668781
42,0	40,0	49	6,0	0,8	35	80623273
42,0	40,0	49	6,0	1,0	43	80623273
42,0	40,0	49	6,0	1,1	47	80623273



Europe

#### Heraeus Deutschland GmbH & Co. KG

Heraeus Performance Products  
 Heraeusstr. 12-14  
 63450 Hanau, Germany  
 Phone +49 6181.35-3740  
 Fax +49 6181.35-8620  
 heraeus-platinum-labware@heraeus.com  
 www.pt-labware.de

United States of America

#### Heraeus Precious Metals

North America Conshohocken LLC  
 Functional Materials  
 24 Union Hill Road  
 W. Conshohocken, PA 12-428-2736, USA  
 Phone +1 484 391 2527  
 Fax +1 610 825 7061  
 heraeus-platinum-labware@heraeus.com  
 www.pt-labware.com