




XRF Platinum Program for Autofluxer & Vulcan




XRF Crucibles suitable for Autofluxer & Vulcan


Diameter internal top (OD) bottom (UD)		Height (H)	Base thick- ness (s)	Weight approx.	Execution	Heraeus drawing
mm	mm	mm	mm	g		
34	20	38	0,34	25	3 pins, with pouring lip	80018374



XRF Crucibles suitable for Autofluxer & Vulcan



Diameter internal top (OD) bottom (UD)		Height (H)	Base thick- ness (s)	Weight approx.	Execution	Heraeus drawing
mm	mm	mm	mm	g		
34	20	38	0,5	43	3 pins on support ring, with pouring lip	80079008
34	20	38	0,5	44	3 pins on support ring with rein- forced rim, with pouring lip	80049875
36	22	38	0,5	51	3 pins on support ring with rein- forced rim, with pouring lip	80623575





XRF Casting dishes suitable for Autofluxer & Vulcan (square version)

Diameter internal		Flange Length (FL)	Height (H)	Base thickness (s)	Weight approx.	Heraeus drawing
top (OD)	bottom (UD)					
mm	mm	mm	mm	mm	g	
31	29	41	3	0,8	27	80623078
31	29	41	3	1,0	34	80623078
34	32	41	3	0,8	31	80623092
34	32	41	3	1,0	38	80623092
36	34	51	3	0,8	45	80623094
36	34	51	3	1,0	58	80623094
41	39	51	3	0,8	45	80623096
41	39	51	3	1,0	58	80623096



XRF Lids suitable for Autofluxer & Vulcan

Diameter	Length (L)	Base thickness (s)	Weight approx.	Execution	Heraeus drawing
mm	mm	mm	g		
45	66	0,3	14	2 pins (straight)	80019941
40	54	0,3	9	2 pins (bent)	81103169



Europe

Heraeus Deutschland GmbH & Co. KG

Heraeus Performance Products
 Heraeusstr. 12-14
 63450 Hanau, Germany
 Phone +49 6181.35-3740
 Fax +49 6181.35-8620
 heraeus-platinum-labware@heraeus.com
 www.pt-labware.de

United States of America

Heraeus Precious Metals

North America Conshohocken LLC
 Functional Materials
 24 Union Hill Road
 W. Conshohocken, PA 12-428-2736, USA
 Phone +1 484 391 2527
 Fax +1 610 825 7061
 heraeus-platinum-labware@heraeus.com
 www.pt-labware.com

